Model CB 170 | CO2 incubators with hot air sterilization and heat sterilizable CO2 sensor

The BINDER CB series CO₂ incubator is the premium class among the CO₂ incubators. It is suitable for all sensitive incubation applications and ensures optimal cell growth. A BINDER CO₂ incubator of the CB series is up to even complex cultivation experiments or individual environments.

BENEFITS

- Safe cell cultures receive maximum protection against contamination
- Reliable constant well-being conditions for the cells
- · Smart simple routine cleaning and convenient operation
- · Economical efficient operation without consumables



Model 170

Model 170

IMPORTANT FEATURES

- Optimized double-pan humidification system with condensation protection
- Tried-and-tested anti-contamination concept with 180°C hot air sterilization
- Hot-air sterilizable CO_2 sensor with single-beam infrared technology
- Seamless inner chamber made of stainless steel with flanges as shelf support
 system
- Fanless interior with Venturi CO2 gas mixing nozzle
- OPTIONAL EQUIPMENT
- O2 control with O2 control ranges: 0.2 to 20 vol.% O2 or 10 to 90 vol.% O2
- Multi-sectional glass doors with special shelves for dividing up the space to accommodate different experiments

- Intuitive touchscreen controller
- Internal data recording, can be read out in open format via USB interface
- Fault diagnosis system with visual and acoustic alarm
- Interfaces: Ethernet, USB, zero-voltage alarm contact
- Fail-safe CO2 system for protecting the pH of cell cultures
- Stacking adapter for secure, space-saving installation
- Gas tank changer for connecting two gas tanks to one or two incubators

| Interior volume [L] | Power supply - unit fuse | Plug* | Version | Model version | ArtNo. |
|------------------------|------------------------------|-----------|-----------------------------|---------------|-----------|
| Model CB 170 | | - | | | |
| | | | Standard | CB170-230V | 9640-0009 |
| | 200230 V 1~ 50/60 Hz -10,0 A | CEE 7/7 | with Access port 30 mm left | CB170-230V-D | 9640-0027 |
| 170 | | | with O2 control 0.2-20% | CB170-230V-0 | 9640-0010 |
| | 100120 V 1~ 50/60 Hz -16,0 A | NEMA 5-20 | Standard | CB170UL-120V | 9640-0011 |

ORDERING INFORMATION

BINDER Best conditions for your success

| Interior volume [L] | Power supply - unit fuse | Plug* | Version | Model version | ArtNo. |
|------------------------|--------------------------|-------|-------------------------------------|----------------|-----------|
| | | | with Access port 30 mm left | CB170UL-120V-D | 9640-0028 |
| | | | with O ₂ control 0.2-20% | CB170UL-120V-O | 9640-0012 |

BINDER Best conditions for your success

| TECHNICAL | DATA |
|-----------|------|
| LCHRICKE | |

| Data | | | | | | |
|--|---|---|---|---|---|--|
| Designation | CB170-230V | CB170-230V-D | CB170-230V-0 | CB170UL-120V | CB170UL-120V-D | CB170UL-120V-0 |
| Article Number | 9640-0009 | 9640-0027 | 9640-0010 | 9640-0011 | 9640-0028 | 9640-0012 |
| Option model | Standard | with Access port 30 mm left | with O2 control 0.2- 20% | Standard | with Access port 30 mm left | with O2 control 0.2- 20% |
| Performance Data Temper | ature | | | | | |
| Temperature range | +4 °C above ambient temperature to 60 °C | +4 °C above ambient temperature to 60 °C | +6 °C above ambient temperature to 60 °C | +4 °C above ambient temperature to 60 °C | +4 °C above ambient temperature to 60 °C | +6 $^{\circ}\text{C}$ above ambient temperature to 60 $^{\circ}\text{C}$ |
| Temperature uniformity at 37°C | 0.3 ± K |
| Temperature fluctuation at 37°C | 0.1 ± K |
| Recovery time after door was opened for 30 s at 37°C | 6 min |
| Performance Data Climate | | | | | | |
| Humidity range | 9095 % RH |
| Gas Data | | | | | | |
| CO2 range | 020 Vol% CO2 |
| CO ₂ measuring technology | IR | IR | IR | IR | IR | IR |
| CO2 recovery time after door was opened for 30 s at 5 vol. % CO2 | 5 min |
| Standard O2 control range | | | 0.220 Vol% O2 | | | 0.220 Vol% O2 |
| O2 control ranges with option: O2 range | | | 1090 Vol% O2 | | | 1090 Vol% O2 |
| O2 recovery time after door was opened for 30 s at 5 vol. % O2 | | | 10 min | | | 10 min |
| Electrical data | | | | | | |
| Rated Voltage | 200230 V | 200230 V | 200230 V | 100120 V | 100120 V | 100120 V |
| Power frequency | 50/60 Hz |
| Nominal power | 1,4 kW |
| Unit fuse | 10,0 A | 10,0 A | 10,0 A | 16,0 A | 16,0 A | 16,0 A |




| Phase (Nominal voltage) | 1- | 1- | 1- | 1~ | 1- | 1~ |
|-----------------------------------|--------------------------------|----------|----------|----------|----------|----------|
| Dimensions and weights | | | | | | |
| Interior volume | 170 L | 170 L | 170 L | 170 L | 170 L | 170 L |
| Net weight of the unit (empty) | 95 kg | 95 kg | 95 kg | 95 kg | 95 kg | 95 kg |
| Load per rack | 10 kg | 10 kg | 10 kg | 10 kg | 10 kg | 10 kg |
| Permitted load | 30 kg | 30 kg | 30 kg | 30 kg | 30 kg | 30 kg |
| Wall clearance back | 100 mm | 100 mm | 100 mm | 100 mm | 100 mm | 100 mm |
| Wall clearance sidewise | 50 mm | 50 mm | 50 mm | 50 mm | 50 mm | 50 mm |
| Housing dimensions not | inc l. fittings and connection | S | | | | |
| Width net | 680 mm | 680 mm | 680 mm | 680 mm | 680 mm | 680 mm |
| Height net | 870 mm | 870 mm | 870 mm | 870 mm | 870 mm | 870 mm |
| Depth net | 715 mm | 715 mm | 715 mm | 715 mm | 715 mm | 715 mm |
| Internal Dimensions | | | | | | |
| Interior width | 560 mm | 560 mm | 560 mm | 560 mm | 560 mm | 560 mm |
| Interior height | 600 mm | 600 mm | 600 mm | 600 mm | 600 mm | 600 mm |
| Interior depth | 500 mm | 500 mm | 500 mm | 500 mm | 500 mm | 500 mm |
| Inner doors | 1 | 1 | 1 | 1 | 1 | 1 |
| Unit doors | 1 | 1 | 1 | 1 | 1 | 1 |
| Environment-specific dat | a | | | | | |
| Sound-pressure level | 41 dB(A) | 41 dB(A) | 41 dB(A) | 41 dB(A) | 41 dB(A) | 41 dB(A) |
| Energy consumption at 37°C | 55 Wh/h | 55 Wh/h | 55 Wh/h | 55 Wh/h | 55 Wh/h | 55 Wh/h |
| Fixtures | | | | | | |
| Number of shelves (std./max.) | 3/6 | 3/6 | 3/6 | 3/6 | 3/6 | 3/6 |

All technical data is specified for unloaded units with standard equipment at an ambient temperature of +22 $^{\circ}C \pm 3 ^{\circ}C$ and a power supply voltage fluctuation of $\pm 10\%$. The temperature data is determined to BINDER factory standard following DIN 12880, observing the recommended wall clearances of 10 % of the height, width, and depth of the inner chamber. All indications are average values, typical for units produced in series. We reserve the right to change technical specifications at any time.



DIMENSIONS Incl. fittings and connections [mm]



OPTIONS

| Designation | Description | CB 170 | * | ArtNo. |
|--|---|--------|----|-----------|
| | back | | | |
| | 30 mm | • | 01 | 8612-0025 |
| Access port with silicone plug | left | | | |
| | 30 mm | • | 01 | 8612-0026 |
| | right | | | |
| | 30 mm | • | 01 | 8612-0027 |
| Analog output 4-20 mA | for temperature and CO_2 values (outputs not adjustable) | • | 02 | 8612-0022 |
| | for temperature, measurement in center of chamber at specified temperature | • | - | 8012-1132 |
| | temperature measurement incl. certificate, 9 measuring points at specified temperature | • | - | 8012-1550 |
| Calibration certificate, temperature | temperature measurement incl. certificate, 15- 18 measuring points at specified temperature | | - | 8012-1571 |
| | temperature measurement incl. certificate and 27 measuring points at specified temperature | • | - | 8012-1592 |
| Calibration certificate, temperature and CO2 | for temperature and CO2, temperature measurement in center of chamber, CO2 measurement performed using test gas at 37 $^\circ$ C and 5 $\%$ CO2 | • | - | 8012-1235 |
| Clibertine continues (2) | for versions with O_2 control; a measuring point in center of chamber at 80 Vol% O_2 | • | - | 8012-2604 |
| Calibration certificate O2 | for versions with O_2 control; a measuring point in center of chamber at 1 Vol% O_2 | • | - | 8012-2603 |
| CELLROLL Set | modular, expandable roller bottle system for cell cultivation, 4 roller bottles | | - | 8612-0023 |
| Door hinged on the left | Unit door and standard glass inner door with door hinged on the left | • | - | 8612-0034 |
| Electric access port | 8-pin, for low voltage with LEMO socket (coverable) and LEMO plug (max. 24 V - 2 A) | | - | 8612-0033 |
| Inner door, divided | Cell therapy compartmentalization, consisting of 6 inner doors and three shelf levels with one partition wall each | | - | 8612-0030 |

Data sheet Model CB 170

BINDER Best conditions for your success

| Designation | Description | CB 170 | * | ArtNo. |
|---------------------|---|--------|---|-----------|
| O2 control range | alternative to the standard range, for controlling O2 concentrations between 10 and 90 vol. %, only in conjunction with option models with O2 control | • | - | 8612-0028 |
| Quick sample access | vewing window (13 x 25 cm) in the inner door for quick access | • | - | 8612-0032 |

ACCESSORIES

| Designation | Description | CB 170 | * | ArtNo. |
|---------------------------------|--|--------|---|-----------|
| | for simple logging and documentation requirements with up to 5 networked units. | | | |
| APT-COM™ 4 BASIC-Edition | version 4, BASIC edition | • | - | 9053-0039 |
| APT-COM™ 4 GLP-Edition | for working under GLP-compliant conditions. Measured values are documented in a tamper-proof way in line with the requirements of FDA Regulation 21 CFR 11. | | | |
| | version 4, GLP edition | • | - | 9053-0042 |
| APT-COM™ 4 PROFESSIONAL-Edition | convenient unit and user management built on the BASIC edition. Suitable for networking up to 100 units. | | | |
| | version 4, PROFESSIONAL edition | • | - | 9053-0040 |
| Base | the base equipped with casters is used for the secure positioning and leveling of a BINDER CO $\mbox{2}$ incubator | | - | 9051-0028 |
| Copper interior fittings | Set consisting of 3 shelves (all-copper) and a water pan (copper-plated) | • | - | 8012-1918 |
| | for CO $_{\rm 2},$ consisting of a gas tank pressure regulator with connection parts and 5-meter hose | • | - | 8012-0014 |
| Gas cylinder connection set | for N_2 , consisting of a gas tank pressure regulator with connection parts and 5-meter hose | • | - | 8012-0016 |
| | for O_2 , consisting of a gas tank pressure regulator with connection parts and 5-meter hose | • | - | 8012-0015 |
| Gas tank changer | external, BINDER Gas Supply Service, for connecting 2 gas tanks (CO ₂ , N ₂ or O ₂), with audible and visual alarms, as well as zero-voltage alarm output | • | - | 8012-2344 |
| pH-neutral detergent | concentrated, for gentle remove of residual contaminants; 1 kg | • | - | 8012-2250 |
| | IQ/OQ/PQ documents - supporting documents for validation performed by customers, according to customer requirements, PQ section added to qualification folder IQ/OQ; parameters: temperature, CO ₂ , O ₂ - or pressure, depending on unit | | | |
| | Digital in PDF format | • | | 7057-0005 |
| Qualification documents | — Hard copy inside folder | • | - | 7007-0005 |

comprehensive unit documentation; parameters: temperature, CO2, O2, pressure, depending on unit



| Designation | Description | CB 170 | * | ArtNo. |
|-----------------------|--|--------|---|-----------|
| | Digital in PDF format | | - | 7057-0001 |
| | Hard copy inside folder | • | - | 7007-0001 |
| Rubber pads | set anti-slip feet | • | - | 8012-0702 |
| Shaker "Professional" | Shaker and controller equipped with LEMO connectors, exclusively for CB with electric access port, 100 - 240 V, 50/60 Hz | • | - | 8012-1772 |
| Shaker "Universal" | Shaker and controller connected via ribbon cable, universally applicable, 100 - 240 V, 50/60 Hz | • | - | 8012-1771 |
| | Stainless steel | • | | 8012-2167 |
| Shelf, perforated | for divided inner door | | | |
| | Stainless steel, divided | • | - | 8012-1888 |
| Stacking adapter | for the thermally isolated stacking of two BINDER CO_2 incubators | • | - | 9051-0035 |
| Stacking frame | stable, vibration-free stacking frame on casters with stop brake for safe stacking of two CB series units | • | - | 9051-0040 |

SERVICES

| Designation | Description | * | ArtNo. |
|--|--|---|-----------|
| Calibration services | | | |
| | Calibration of one (1) test temperature specified by the user in center of chamber, including certificate | - | DL30-0140 |
| Calibration certificate, temperature | Extension of calibration of one (1) additional test temperature specified by the user in the center of the usable space, including certificate | _ | DL30-0102 |
| CO2 calibration | CO2 calibration at a specified value, measurements are taken with analyzed test gas at 5%, including certificate | _ | DL30-0401 |
| O2 calibration | O_2 calibration with 1 measuring point in center of chamber at a specified O_2 percentage value, including certificate | _ | DL30-0402 |
| Temperature measurement, 9 measuring points | Temperature measurement with 9 measuring points with a set value specified by the user, including certificate | _ | DL30-0109 |
| Temperature measurement, 18 measuring points | Temperature measurement with 18 measuring points with a set value specified by the user, including certificate | - | DL30-0118 |
| Temperature measurement, 27 measuring points | Temperature measurement with 27 measuring points with a set value specified by the user, including certificate | - | DL30-0127 |
| Installation services | | | |
| Unit commissioning | Connect the unit to the customer-side connections (electricity, water, wastewater, gas), basic functions check, brief operating instructions. (excl.: unpacking, setup, controller instructions, programming, installation work) | - | DL10-0110 |
| Unit instructions | Instruction regarding operating principle and basic functions of the unit, operation of the control electronics including programming | - | DL10-0510 |
| Maintenance contracts | | | |

Data sheet Model CB 170

BINDER Best conditions for your success

| escription | * | ArtNo. |
|---|---|---|
| intenance service as contractually arread visual inspection of mechanical and | | |
| ectrical components, check of control response, 20% discount on spare parts | - | DL20-0710 |
| intenance service as contractually agreed, visual inspection of mechanical and ectrical components, check of control response, 20% discount on spare parts, sting of all key functions, replacement of wear parts, calibration of one test mperature specified by the user in the center of the usable space, including rtificate | - | DL20-0975 |
| intenance service as contractually agreed, visual inspection of mechanical and ectrical components, check of control response, 20% discount on spare parts, sting of all key functions, calibration of one test temperature specified by the er in the center of the usable space, without certificate | _ | DL20-0875 |
| | | |
| ne-off maintenance service in accordance with maintenance schedule. Visual spection of mechanical and electrical components, testing of all key functions. libration of a test temperature specified by the user in center of usable space thout certificate | - | DL20-0608 |
| | | |
| ecution of IQ/OQ in accordance with qualification folder | - | DL41-0200 |
| ecution of IQ/OQ/PQ in accordance with qualification folder | - | DL44-0500 |
| | | |
| e warranty is extended by 1 year from the delivery date, wear parts are excluded | - | DL50-0020 |
| | intenance service as contractually agreed, visual inspection of mechanical and ctrical components, check of control response, 20% discount on spare parts, ting of all key functions, replacement of wear parts, calibration of one test inperature specified by the user in the center of the usable space, including rtificate intenance service as contractually agreed, visual inspection of mechanical and ctrical components, check of control response, 20% discount on spare parts, ting of all key functions, calibration of one test temperature specified by the er in the center of the usable space, without certificate me-off maintenance service in accordance with maintenance schedule. Visual pection of mechanical and electrical components, testing of all key functions. Libration of a test temperature specified by the user in center of usable space thout certificate | intenance service as contractually agreed, visual inspection of mechanical and ctrical components, check of control response, 20% discount on spare parts, ting of all key functions, replacement of wear parts, calibration of one test - nperature specified by the user in the center of the usable space, including rtificate |

NOTES

- 01 Condensation may occur in the area around the access port. Access ports may be placed in custom locations for an additional charge.
- 02 UL mark is not granted when this option is used.

BINDER GmbH Tuttlingen, Germany TEL +49 7462 2005 0 info@binder-world.com www.binder-world.com

BINDER Asia Pacific (Hong Kong) Ltd. Kowloon, Hong Kong, P.R. China TEL +852 39070500 asia@binder-world.com www.binder-world.com BINDER Inc. Bohemia, NY, USA TEL +1 631 224 4340 usa@binder-world.com www.binder-world.us BINDER Environmental Testing Equipment (Shanghai) Co., Ltd. Shanghai, P.R. China TEL +86 21 685 808 25 china@binder-world.com www.binder-world.com