



## Model MKT 240 | Dynamic climate chambers for rapid temperature changes with extended low temperature range

Temperature ranges between -70 °C and 180 °C, with the added benefit of natural simulation, are what make the BINDER MKT series so unique. At the same time, this environmental simulation chamber fulfills the very highest precision and performance requirements for cyclical temperature tests.

#### **BENEFITS**

- Homogeneous climate conditions thanks to APT.line<sup>™</sup> technology
- Comprehensive programming and data acquisition
- · Large heated viewing window



Model 240



Model 240

#### IMPORTANT FEATURES

- Temperature range: -70 °C to +180 °C
- 4 zero-voltage relay contacts
- APT.line<sup>™</sup> preheating chamber technology
- Programmable condensation protection for test material
- Heated viewing window with LED interior lighting
- BINDER Multi Management Software APT-COM™ Basic Edition
- · Troubleshooting system with visual and audible alarms
- · Intuitive touchscreen controller with time-segment and real-time programming · Inner chamber made of stainless steel
- Internal data logger, measured values can be read out in open format via USB CFC-free refrigerants R-452A and R-23
- Access port with silicone plug: 50 mm, left
- Class 2 independent adjustable temperature safety device (DIN 12880) with visual alarm

- 4 castors, two with brakes
- Computer interface: Ethernet
- Adjustable ramp function
- Integrated chart recorder
- Real-time clock
- Door heating
- 1 stainless steel rack

- Cooling with cascade compressor cooling unit

#### ORDERING INFORMATION

Interior volume [L]	Power supply - unit fuse	Plug*	Version	Model version	ArtNo.
Model MKT 240					
220	400 V 3~/N/PF 50 H7 - 16 O A	CEE 16 / 6 H plug 5- pin	Standard	MKT240-400V	9020-0386
228	480 V 3~ 60 Hz -16,0 A	CEE 16 / 6 H plug 5- pin	with voltage and frequency converter	MKT240-480V-C	9020-0364



#### TECHNICAL DATA


TECHNICAL DATA		
Data		
Designation	MKT240-400V	MKT240-480V-C
Article Number	9020-0386	9020-0364
Option model	Standard	with voltage and frequency converter
Performance Data Temperature		
Temperature range	-70180 ℃	-70180 ℃
Temperature uniformity dependent on set value	0.1. 1.0 ± K	0.1. 1.0 ± K
Temperature fluctuation dependent on set value	0.1. 0.4 ± K	0.1. 0.4 ± K
Average heating-up rate according to IEC 60068-3-5	5.0 K/min	5.0 K/min
Cooling down time from 180°C to -70°C	95 min	95 min
Average cooling down time according to IEC 60068-3-5	4.2 K/min	4.2 K/min
Electrical data		
Rated Voltage	400 V	480 V
Power frequency	50 Hz	60 Hz
Nominal power	6,5 kW	6,5 kW
Unit fuse	16,0 A	16,0 A
Phase (Nominal voltage)	3~/N/PE	3~
Dimensions and weights		
Interior volume	228 L	228 L
Net weight of the unit (empty)	380 kg	380 kg
Load per rack	30 kg	30 kg
Permitted load	70 kg	70 kg
Wall clearance back	300 mm	300 mm
Wall clearance sidewise	200 mm	200 mm
Housing dimensions not incl. fittings and connections		
Width net	1,115 mm	1,115 mm
Height net	1,935 mm	1,935 mm
Depth net	925 mm	925 mm
Internal Dimensions		
Interior width	735 mm	735 mm




Interior height	700 mm	700 mm			
Interior depth	443 mm	443 mm			
Unit doors	1	1			
Environment-specific data					
Sound-pressure level	64 dB(A)	67 dB(A)			
Energy consumption at 20°C	500 Wh/h	500 Wh/h			
Fixtures					
Number of shelves (std./max.)	1/6	1/6			


All technical data is specified for unloaded units with standard equipment at an ambient temperature of  $\pm 2^{\circ}$  C  $\pm 3^{\circ}$  C and a power supply voltage fluctuation of  $\pm 10^{\circ}$  %. The temperature data is determined in accordance to BINDER factory standard following DIN 12880, observing the recommended wall clearances of 10 % of the height, width, and depth of the inner chamber. Technical data refers to 100 % fan speed. All indications are average values, typical for units produced in series. We reserve the right to change technical specifications at any time.

## DIMENSIONS Incl. fittings and connections [mm]










Heating up rate

Cooling down rate



Heat compensation with a set value of -70°C  $\,$ 

Heat compensation chart

## **OPTIONS**

Designation	Description	MKT 240	*	ArtNo.
	left			
	30 mm	•	01	8012-1324
	50 mm	•	01	8012-1330
Access part with cilicons plug	80 mm	•	01	8012-1336
Access port with silicone plug	100 mm	•	01	8012-1544
	125 mm	• 01	8012-1353	
	right			
	30 mm	•	01	8012-1321



Designation	Description	MKT 240	*	ArtNo.
	50 mm	•	01	8012-1327
	80 mm	•	01	8012-1333
	100 mm	•	01	8012-1541
	125 mm	•	01	8012-1350
	top			
	80 mm	•	01	8012-1538
	100 mm	•	01	8012-1532
	125 mm	•	01	8012-1535
Analog output 4-20 mA	for temperature values (output not adjustable)	•	-	8012-1785
Calibration certificate, expanded	for temperature; for extending the measurement in center of chamber to include another test temperature	•	-	8012-1125
	for temperature, measurement in center of chamber at specified temperature	•	-	8012-1144
Calibration certificate,	temperature measurement incl. certificate, 9 measuring points at specified temperature	•	-	8012-1561
temperature	temperature measurement incl. certificate, 15- 18 measuring points at specified temperature	•	-	8012-1582
	temperature measurement incl. certificate and 27 measuring points at specified temperature	•	-	8012-1603
Compressed air connection	for the connection to an existing pressurized air network	•	-	8012-1788
Door lock	lockable door handle	•	-	8012-1784
Dry-air purge	uncontrolled, incl. connection; for preventing condensation	•	-	8012-1867
Notch-type access port	notch-type access port in door, 100 x 35 mm	•	-	8012-1855
Pt 100 temperature sensor	additional flexible Pt 100, interior, for displaying the temperature on the unit display	•	-	8012-1790
RS 485 interface, 2-wire	Additional serial interface can be used parallel to Ethernet, for Multi Management Software APT-COM™	•	-	8012-1770
Temperature safety device class 2	with visual alarm (DIN 12880)	•	-	8012-1794

### **ACCESSORIES**

Designation	Description	MKT 240	*	ArtNo.
APT-COM™ 4 GLP-Edition	for working under GLP-compliant conditions. Measured values are documented in a tamper-proof way in line with the requirements of FDA Regulation 21 CFR $11$ .			
	version 4, GLP edition	•	_	9053-0042



Designation	Description		MKT 240		*	ArtNo.
APT-COM™ 4 PROFESSIONAL- to 100 units. Edition		nd user management built on the BASIC edition. Suitable for networking up				
		ESSIONAL edition	•		-	9053-0040
pH-neutral detergent	concentrated, for	gentle remove of residual contaminants; 1 kg	•		-	8012-2250
	according to cust	nents – supporting documents for validation performed by customers, omer requirements, PQ section added to qualification folder IQ/OQ; operature, CO <sub>2</sub> , O <sub>2</sub> – or pressure, depending on unit				
	Digital in PDF fo	ormat _	•		-	7057-0005
Ovelifienties de suscepte	Hard copy insid	e folder	•		-	7007-0005
Qualification documents	of: IQ/OQ check	s – supporting documents for validation performed by customers, consisting lists incl. calibration guide and comprehensive unit documentation; perature, CO <sub>2</sub> , O <sub>2</sub> , pressure, depending on unit				
	Digital in PDF fo	ormat	•		-	7057-0001
	Hard copy insid	e folder	•		-	7007-0001
Rack	stainless steel		•		-	8012-2053
Rack, reinforced	stainless steel, with fasteners (1 set of 4)		•		-	8012-0605
Rack accessories	fasteners (1 set of 4) for additional security of racks		•		-	8012-2280
	RS 422 cable set distributor	and RS 485 / RS 422 interface converter for connection to 10-way plug				
RS 485 / RS 422 interface converter	115 V option mo	del	•		-	8012-0599
	230 V option m	odel –	•		-	8012-0589
Shelf, perforated	Stainless steel		•		-	8012-2251
SERVICES						
Designation		Description		*	Art	No.
Calibration services						
Calibration certificate, temperature		Calibration of one (1) test temperature specified by the user in center of chamber including certificate	er,	-	DL	30-0130
		Extension of calibration of one (1) additional test temperature specified by the user in the center of the usable space, including certificate		-	DL	30-0102
Temperature measurement, 9 measuring points  Temperature measurement with 9 measuring points with a set value specifies the user, including certificate		ру		DL	30-0109	
Temperature measurement, 18 measuring points the user, including certificate		Temperature measurement with 18 measuring points with a set value specified the user, including certificate	by	-	DL3	30-0118
Temperature measurement, 27 measuring points  Temperature measurement with 27 measuring points with a set value specific the user, including certificate		by	_	DLS	30-0127	



Designation	Description	*	ArtNo.
Installation services			
Unit commissioning	Connect the unit to the customer-side connections (electricity, water, wastewater, gas), basic functions check, brief operating instructions. (excl.: unpacking, setup, controller instructions, programming, installation work)	-	DL10-0300
Unit instructions	Instruction regarding operating principle and basic functions of the unit, operation of the control electronics including programming	-	DL10-0700
Maintenance contracts			
BRONZE 3-year maintenance contract	Maintenance service as contractually agreed, visual inspection of mechanical and electrical components, check of control response, 20% discount on spare parts	_	DL20-0710
GOLD 3-year maintenance contract	Maintenance service as contractually agreed, visual inspection of mechanical and electrical components, check of control response, 20% discount on spare parts, testing of all key functions, replacement of wear parts, calibration of one temperature/humidity/pressure value, including certificate	-	DL20-0930
SILVER 3-year maintenance contract	Maintenance service as contractually agreed, visual inspection of mechanical and electrical components, check of control response, 20% discount on spare parts, testing of all key functions, calibration of one test temperature specified by the user in the center of the usable space, without certificate	-	DL20-0830
Maintenance services			
Leakage test	Testing and logging of the leak-tightness of systems containing hydrofluorocarbons in accordance with Regulation EC842/2006, incl. test book	-	DL00-0034
Maintenance	One-off maintenance service in accordance with maintenance schedule. Visual inspection of mechanical and electrical components, testing of all key functions. Calibration of a test temperature specified by the user in center of usable space without certificate	-	DL20-0607
Validation services			
Execution of IQ/OQ	Execution of IQ/OQ in accordance with qualification folder	-	DL41-0200
Execution of IQ/OQ/PQ	Execution of IQ/OQ/PQ in accordance with qualification folder	_	DL44-0500
Warranty service			
1-year warranty extension	The warranty is extended by 1 year from the delivery date, wear parts are excluded	-	DL50-0040

# BINDER Best conditions for your success

#### **NOTES**

On Condensation may occur in the area around the access port. Access ports may be placed in custom locations for an additional charge.

BINDER GmbH Tuttlingen, Germany TEL +49 7462 2005 0 info@binder-world.com www.binder-world.com BINDER Inc. Bohemia, NY, USA TEL +1 631 224 4340 usa@binder-world.com www.binder-world.us BINDER Environmental Testing Equipment (Shanghai) Co., Ltd. Shanghai, P.R. China TEL +86 21 685 808 25 china@binder-world.com www.binder-world.com

BINDER Asia Pacific (Hong Kong) Ltd. Kowloon, Hong Kong, P.R. China TEL +852 39070500 asia@binder-world.com www.binder-world.com