

Klik hier voor meer informatie op vonmarcken.nl

# PRODUCT

## -60°C LOW TEMP. FREEZERS

The VT low temperature freezers provides the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning makes the VT LT freezers an ideal solution for laboratories and hospitals.

Nature R refrigerant, cyclopentane insulation combined with recyclable materials, ensures that the freezers are extremely eco-friendly.

### SPECIFICATIONS

Gross volume, litres (cu. ft.)	296 (10.5)
Internal storage capacity, litres (cu. ft.)	284 (10.0)
Temperature range (16°C-30°C AMB)	-30°C to -60°C
Temperature probes	1
Compressor	Embraco
Refrigerant	Nature R
Noise level, dB	55
Power consumption, kWh / 24h	6,2
Power consumption, kWh / year	2.263
Voltage, V	220-240
Frequency, Hz	50



### FEATURES

Castors	YES
No. of baskets	3
Defrost	Manual
Battery Backup	NO
USB	NO
Data connection	NO/NC
Digital controller	XR30CX
Data logging, # logs	NO
Analogue chart recorder	NO
Max no. of standard vials, pcs.	17.496

### DIMENSIONS

Height, mm (inches)	850 (33.5)
Height w. castors, mm (inches)	890 (35.1)
Width, mm (inches)	1.260 (28.3)
Depth, mm (inches)	695 (27.4)
Gross weight, kg (lbs)	82 (181)
Net weight, kg (lbs)	76 (168)

### LOADING QUANTITIES

Qty. per 20' / 40' container	24 / 54
Qty. per truck (EU standard)	60



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.



- · State-of-the-art unit designed and manufactured by A/S Vestfrost in Denmark
- Unit optimized for optimal footprint

- · Eco-friendly, HCFC-free Nature R refrigerant
- · Efficient cyclopentane insulation
- · Low energy consumption

- · Digital XR30CX controller with display
- · Prepared for external logging equipment
- PT1000 internal temperature measuring probe
- · Top-quality castors for explicit movability
- · Integrated lock to prevent unauthorized access
- · High/low temperature alarm

- · User-friendly and maintenance-free
- · Single compressor system
- · Suitable for rack systems
- · ISO 9001, 14001 & 18001 OHSAS Certified



-60°C LOW TEMP. FREEZERS

**VT 307** 



PRECISION digital XR30CX controller with display, which is pre-



 $\ensuremath{\text{TOP-QUALITY}}$  castors eases repositioning of the VT line units (Optional on some VT units).



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.





**BIOMEDICAL SOLUTIONS** 

# **DIMENSIONS**

# -60°C LOW TEMP. FREEZER

250





### DIMENSIONS

Height, mm (inches)	850 (33,46)
Height w. castors, mm (inches)	890 (35,04)
Width, mm (inches)	1.260 (49,61)
Depth, mm (inches)	600 (23,62)
Gross weight, kg (lbs)	82 (181)

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

All measurements in mm

# DIMENSIONS

## -60°C LOW TEMP. FREEZER



### **INTERNAL DIMENSIONS**

Height, mm (inches)	635 (25)
Width, mm (inches)	1.100 (43,31)
Depth, mm (inches)	440 (17,32)
Insulation, mm	80
Insulation type	PUR
Innerliner material	Coated steel
Max. AMB humidity, %	70
Internal humidity, %	<5

## LOADING CAPACITY

Max loading, kg (lbs)	200 (440)
Max no. of standard vials, pcs.	17.496
Max loading in baskets, kg (lbs)	168 (370)

## **UNIT FEATURES**


Porthole for external probe, mm $\ensuremath{\mathcal{Q}}$	12
Castors (with brakes)	4 (0)
Inner lids	0
Lid vacuum release	N/A

### **BASKET FEATURES**

Basket material	Coated steel
Basket type	Medium
Quantity	3
Maximum basket load, kg (lbs)	14,50 (31,96)
Weight, kg (lbs)	0,390 (0,86)



A/S Vestfrost Falkevej 12 DK-6705 Esbjerg Ø www.vestfrostsolutions.com



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

# RACKS

## VT 307 -60°C LOW TEMP. FREEZERS

### STORAGE CAPACITY: RACKS FOR 50 MM BOXES

VF Product no.	No. of boxes	No. of racks	Total no. of boxes	
VF311200	6	3	18	
VF311220	11	18	198	
Total, ex dividers			216	- AND

## STORAGE CAPACITY: RACKS FOR 50 MM BOXES (INCL. DIVIDERS)

VF Product no.	No. of boxes	No. of racks	Total no. of boxes
VF312200	6	3	18
VF312220	11	18	198
Total, vials			17.496

### STORAGE CAPACITY: RACKS FOR 75 MM BOXES

VF Product no.	No. of boxes	No. of racks	Total no. of boxes	at to the
VF311300	4	3	12	
VF311320	8	18	144	
Total, ex dividers			156	

### STORAGE CAPACITY: RACKS FOR 75 MM BOXES (INCL. DIVIDERS)

VF Product no.	No. of boxes	No. of racks	Total no. of boxes
VF312300	4	3	12
VF312320	8	18	144
Total, vials			12.636

### STORAGE CAPACITY: RACKS FOR 100 MM BOXES

VF Product no.	No. of boxes	No. of racks	Total no. of boxes	
VF311400	3	3	9	
VF311420	5	18	90	
Total, ex dividers			99	

STORAGE CAPACITY: RACKS FOR 100 MM BOXES (INCL. DIVIDERS)

VF Product no.	No. of boxes	No. of racks	Total no. of boxes
VF312400	3	3	9
VF312420	5	18	90
Total, vials			8.019



